




Combining Rapid Isomer Separations and Physicochemical Property Predictions for Drug Molecules with Differential Mobility Spectrometry

Chang Liu,¹ J. C. Yves Le Blanc,¹ Jefry Shields,² Hui Zhang,² John S. Janiszewski,² Christian Ieritano,³ Luke Melo,³ Evan Shepherson,³ Mitch Verbuyst,³ Moaraj Hasan,³ Dalia Naser,³ W. Scott Hopkins,³ J. Larry Campbell,¹ Tim L. Hoffman¹


¹SCIEX, Concord, ON; ²Pfizer, Groton, CT; ³University of Waterloo, Waterloo, ON

ASMS 2016, San Antonio, June 6th, 2016

Drug designers are keen to speed up drug development



Good cell permeability




Good drug transport



Good drug candidate


Blausen.com staff. "Blausen gallery 2014". Wikiversity Journal of Medicine.

Differential mobility spectrometry (DMS)



**SCIEX SelexION® Differential Mobility Spectrometer with
QTRAP® 5500 Mass Spectrometer System**


Differential mobility spectrometry (DMS)



The magnitude of the compensation voltage is a measure of the ion-solvent interaction.

Three types of DMS behaviors



Stronger ion-solvent interaction




Isomer separation with DMS

Steric effects can be a dominating factor in isomer separation

Volume 140 | Number 20 | 21 October 2015 | Pages 6761–7032




No modifier (N₂ Only)



Type C


Water



Type B

Analyst

www.rsc.org/analyst



Themed issue: Ion Mobility Mass Spectrometry

ISSN 0003-2654






PAPER
W. Scott Hopkins, J. Larry Campbell et al.
Using differential mobility spectrometry to measure ion solvation: an examination of the roles of solvents and ionic structures in separating quinoline-based drugs




Liu et al., *Analyst*, **2015**, 140, 6897-6903.

Isomer separation with DMS

Steric effects can be a dominating factor in isomer separation in the presence of solvent (e.g., water, methanol)




Water
Type B




Isomer separation with DMS: Next level drug molecules – Quinolin-8-ols


5-substituted
2-methyl-quinolin-8-ol



6-substituted
2-methyl-quinolin-8-ol




7-substituted
2-methyl-quinolin-8-ol





X = -OCH₃, -CH₃, -H, -Br, -Cl, -F, -NO₂, or -CN



DMS ionograms reveal ion/solvent interaction strength: 5-, 6-, and 7-substituted 2-methyl-8-hydroxyquinolines and MeOH





5-substituted
2-methyl-quinolin-8-ol



6-substituted
2-methyl-quinolin-8-ol




7-substituted
2-methyl-quinolin-8-ol





Brown and Okamoto

J. Am. Chem. Soc. 1958, 80, 4979


DMS correlation with calculated binding energy



5-substituted
2-methyl-quinolin-8-ol







6-substituted
2-methyl-quinolin-8-ol



Wales and Doye, *J. Phys. Chem. A* 1997, 101, 5111; Hopkins et al., *J. Phys. Chem. A* 2013, 117, 10714

Cell permeability (RRCK) data among isomeric drug molecules also agrees with DMS



Conclusions

- DMS is a powerful platform for isomer separations, for studying the ion-solvent interaction strength, and ...
- DMS shows promise as a fast approach to predict the cell permeability (RRCK) of drug candidates.
- Data strongly correlates to
 - calculated binding energies with solvent molecules
 - physical organic constants from the 1950s
 - RRCK data among isomeric drug molecules
 - pKa
- Both steric effect and electronic effect influence the ion-solvent interaction.
 - Steric effect: weak ion-solvent interaction when a bulk substituent is close to the charge site
 - Electronic effect: resonance structures disperse charge density and reduce ion-solvent interactions

Acknowledgements

- Drs. Bradley Schneider and Thomas Covey (SCIEX)
- High performance computing support from the SHARCNET consortium of Compute Canada
- Prof. Terry McMahon (University of Waterloo)
- NSERC (ENGAGE & ENGAGE Plus grants)
- RSC for the permission to reuse published figures.

Trademarks/Licensing

AB Sciex is doing business as SCIEX.

For Research Use Only. Not for use in diagnostic procedures.

© 2016 AB Sciex. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license.

Document number: RUO-MKT11-4152.