

Isolation of Isomeric Drug Metabolites by Differential Mobility Spectrometry: A Proof of Concept Study using Caffeine

Christian Ieritano, Bradley Howie, Larry J. Campbell, W. Scott Hopkins 1*

¹University of Waterloo, 200 University Ave W., Waterloo, ON, N2L 3G1, ^bSCIEX, 71 Four Vally Drive, Concord, ON, L4K 4V8

Abstract

Caffeine is one of the most widely consumed psychoactive established pharmacokinetic substances Differential Mobility Spectrometry- Mass Spectrometry (DMS-MS) has been used in the past to separate complicated mixtures of analytes. We envision a proof of concept study whereby DMS can be used to isolate metabolites of caffeine and their tautomer's to quantify their relative abundance. This could prove useful in determining physicochemical properties of metabolites in clinical trials on new drug species.




Figure 1: Metabolic pathway of

Experimental Methods

DMS was performed on a SELEXion equipped with a 5500 QTRAP.

- Compensation voltages (CV) were ramped at particular separation voltage (SV) to generate an ionogram
- Collision induced dissociation (CID) used to generate mass spectra of breakdown products
- Dispersion plots were created by taking the maximum CV of each peak and plotting it against SV
- H_2O modifier was infused into the N_2 carrier gas at 1.5% (v/v)




Figure 3: Three behaviours observed in the DMS cell²

Computational Methods

Computational portion was completed using Gaussian 09 Caffeine and its metabolites were optimized at the B3LYP/6-311++g(d,p) level of theory Each possible protonation site belonging to each tautomeric form were explored Solvent binding energies were determined employing counterpoise and empirical dispersion corrected






Figure 4: Global minima of tautomeric [Paraxanthine] and [Paraxanthine·H₂O]

Experimental Results



- Breakdown curves were constructed which indicate unique tautomeric protonation sites (data not shown)
- Primary three Metabolites in Dry N₂ exhibited Type C and Type B clustering behavior
- Primary three Metabolites in H₂O exhibited Type B and A clustering behavior
- In H₂O modifier it appears that many tautomeric forms exist and this will further be investigated

Computational Results



- The three bare protonated paraxanthine tautomer's that were DFT-optimized ranked in shown figure 6 based on their relative energies
- Protonated paraxanthine tautomer's with H₂O are shown ranked based on their relative energies
- the relative of solvent binding energies found from the protonated tautomer's complexed with H₂O are shown ranked

Figure 6: i) Relative energies of paraxanthine-keto and enol ii) Relative energies of paraxanthine-keto and enol in H_2O iii) H_2O solvent binding energies

Experimental Directions

- DMS traces and calculations of the remaining metabolites will be obtained
- To explore the validity of DMS to isolate caffeine metabolites, technique utilizing coated blade spray (CBS) will be employed




Figure 7: CBS-DMS-MS schematic

References and Acknowledgments

(1) Rybak, M.E, Pao, C., Anal. Bioanl. Chem., **2014**, 406, 771-784. (2) Hopkins, W.S, Mol. Phys., 2015, 113 (21),

3151-3158.









UNIVERSITY OF WATERLOO **FACULTY OF SCIENCE**



